MECHANICS / OSCILLATIONS

UE1050550

POHL'S TORSION PENDULUM



EXPERIMENT PROCEDURE

- Measure the amplitude of forced oscillations as a function of the excitation frequency for various degrees of damping.
- Observe the phase shift between the excitation and the actual oscillation for excitation frequencies which are very small and other which are very large.

OBJECTIVE Measurement and analysis of forced harmonic rotary oscillation

SUMMARY

0

Pohl's wheel or rotating (torsional) pendulum allows for the investigation of forced harmonic rotary oscillation. For this purpose, the oscillating system is connected to an excitation linkage which is driven by an adjustable-speed DC motor so that the restoring spring periodically extends and compresses. In this experiment the amplitude is measured as a function of the excitation frequency for various degrees of damping and the phase shift between the excitation and the actual oscillation is observed.

REQUIRED APPARATUS

uantity	Description	Number	
1	Pohl's Torsion Pendulum	1002956	
1	Mechanical Stopwatch, 15 min	1003369	
1	Plug In Power Supply 24 V, 0.7 A (230 V, 50/60 Hz)	1000681 or	
	Plug In Power Supply 24 V, 0.7 A (115V, 50/60 Hz)	1000680	
1	DC Power Supply 0 – 20 V, 0 – 5 A (230 V, 50/60 Hz)	1003312 or	
	DC Power Supply 0 – 20 V, 0 – 5 A (115 V, 50/60 Hz)	1003311	
2	Analogue Multimeter AM50	1003073	
1	Set of 15 Safety Experiment Leads, 75 cm	1002843	

BASIC PRINCIPLES

Pohl's wheel or rotating (torsional) pendulum allows for the investigation of forced harmonic rotary oscillation. For this purpose, the oscillating system is connected to an excitation linkage which is driven by an adjustable-speed DC motor so that the restoring spring periodically extends and compresses.

The equation of motion for this system is as follows

(1)
$$\frac{d^2\varphi}{dt^2} + 2\cdot\delta\cdot\frac{d\varphi}{dt} + \omega_0^2\cdot\varphi = A\cdot\cos(\omega_{\rm E}\cdot t)$$

where

J: moment of inertia D: spring constant k: damping coefficient M_0 : amplitude of external torque $\omega_{\rm F}$: angular frequency of external torque

 $\delta = \frac{k}{2I}, \ \omega_0^2 = \frac{D}{I}, \ A = \frac{M_0}{I}$

The solution to this equation is composed of a uniform and a non-uniform component. The uniform component is equivalent to damped simple harmonic motion, as investigated in experiment UE1050500. This decreases exponentially over time and can be neglected by comparison with the nonuniform component after a short period of settling. The non-uniform component

(2)
$$\varphi(t) = \varphi_{\rm E} \cdot \cos(\omega_{\rm E} \cdot t - \psi_{\rm E})$$

is linked to the external torque, however, and remains non-negligible as long as that torque is present: Its amplitude is as follows:

(3)
$$\varphi_{\rm E} = \frac{A_0}{\sqrt{\left(\omega_0^2 - \omega_{\rm E}^2\right)^2 + 4 \cdot \delta^2 \cdot \omega_{\rm E}^2}}$$

This becomes increasingly high the closer the excitation frequency $\omega_{\rm F}$ is to the intrinsic resonant frequency ω_0 of the rotating pendulum. Resonance is said to occur when $\omega_{\rm F} = \omega_0$.

The phase shift is shown below:

(4)
$$\psi_{\rm E} = \arctan\left(\frac{2\cdot\delta\cdot\omega_{\rm E}}{\omega_0^2-\omega_{\rm E}^2}\right)$$

This indicates that the deflection of the pendulum lags behind the excitation. For low frequencies it is close to zero but as the frequency increases, it rises, reaching 90° at the resonant frequency. For very high excitation frequencies, the excitation and oscillation frequencies end up being 180° out of phase.



EVALUATION

The amplitudes of the damped oscillations are plotted against the excitation frequency. This results in a selection of curves which can be described by equation (4) as long as the appropriate damping parameter δ is chosen.

There will be slight deviations from the damping values measured in experiment UE1050500. This is mainly due to the fact that the force of friction is not exactly proportional to the speed as assumed here.



Fig. 1: Resonance curves for various degrees of damping