



SAFE EXPERIMENTING IN CLASS MODERN MEASURING INSTRUMENTS

3bscientific.com

Power supplies • Multimeters • Function generators • Oscilloscopes • Digital thermometers • Sound level meters • pH meters •
Geiger counters • Measurement amplifiers



Dear Customers,

Carrying out experiments on one's own is a fascinating and important part of any science lesson. Of course, though, it is essential never to lose sight of students' safety. Student experiments on open electric circuits, for instance, are particularly demanding in terms of electrical safety.

In this brochure we have put together a selection of electricity power supplies for you, each of which is equipped with a safety transformer compliant with European standard EN 61558-2-6 and guarantees safe isolation between the mains supply and output circuitry.

It goes without saying that our measuring instruments also meet the utmost demands for electrical safety. For analogue voltage and current measurements, we recommend using our ESCOLA range of multimeters and our new demo-multimeter or, for measuring very small charges, currents and voltages, our electrometer or the new U-series measurement amplifier.

Moreover, we also offer you a selection of modern digital measuring instruments. Among the new items in the range are two inexpensive digital oscilloscopes, a thermal imaging camera, an infra-red thermometer, two pH meters and a laser range finder. See for yourselves the many and varied uses to which they can be put.

Our whole project range for physics can be found in our 3B Scientific Physics catalogue or online at 3bscientific.com.

If you do not have the catalogue to hand at the moment, that is no problem. We will be happy to send you your own personal example. Alternatively you can simply visit us online. All of our catalogues are available for download in PDF format.

We hope you have a lot of fun in discovering our selection of new products. We would always be happy to hear from you about your needs or suggestions and look forward to your ordering from us. Our expert team is of course at your disposal at any time for personal advice.

Yours faithfully,

Dr. Johannes Selbach
Head of Product Management Natural Sciences



Delivery free for inland deliveries
of value in excess of 200 €!

*10% discount for all products in
this brochure, promotion closes
31/07/17



Your 10% discount is this easy to get:

1. Visit us at 3bscientific.com
2. Select your product and put it in your shopping basket.
3. Enter the promotion code **PF71** when you are ordering the shopping basket or state it along with your telephone order at +44 1934 425333.

CONTENTS

- 3 Power supplies
- 6 Multimeters
- 8 Function generators
- 9 Digital oscilloscopes
- 10 Thermal imaging camera / Digital thermometers
- 12 Sound level meters
- 13 pH meters / Miscellaneous
- 14 Geiger counters
- 15 Measurement amplifiers

> Ideal for electron tubes!



1013412

> Ideal for experiments on electrostatics!



1019234

> **NEW!**

High-Voltage Power Supply 10 kV (230 V, 50/60 Hz)

Universally applicable, floating high-voltage source for experiments on electrostatics or for operating electron tubes. With built-in transformer resistance to external voltage to provide the heating voltage for electron tubes. Continuously adjustable high-voltage source, which is safe to touch, with passive current limiting and digital voltage display. A temperature-controlled fan protects the equipment from overheating.

- Current limiting to protect against contact with high voltages
- Safety transformer conforming to EN 61558-2-6
- Safe isolation between power supply and output circuits

High-voltage output: 0 – 10000 V DC, max. 2 mA, floating
 Heater voltage output: 6.3 V AC, max. 3 A, resistance to voltage up to 10 kV

1019234 753,30€

More Power Supplies at 3bscientific.com!

DC Power Supply 450 V (230 V, 50/60 Hz)

Power supply with three outputs for supplying electrical power in experiments using an electrometer (1001025).

- Current limiting to protect against contact with high voltages
- Safety transformer conforming to EN 61558-2-6
- Safe isolation between power supply and output circuits

Output 1:

Voltage: 0 – 450 V DC
 Max. current: 10 μ A

Output 2:

Voltage: 1.2 – 12 V DC
 Max. current: 100 mA

Output 3:

Voltage: 0 – 12 V AC
 Max. current: 10 mA
 Dimensions: 250x100x160 mm³ approx.

Weight: 0.8 kg approx.

1008535 265,50€

ADVANTAGES

- + Very high-quality and extremely lightweight equipment in modern casing
- + 3-digit digital display for high voltage
- + Regulated high-voltage output not dependent on mains voltage
- + High-voltage-proof supply for heater voltage
- + No need to change fuses

High-Voltage Power Supply E (230 V, 50/60 Hz)

Universally applicable, floating high-voltage source for electrostatic experiments and for operating spectral tubes, gas discharge tubes and electron tubes. With built-in transformer resistance to external voltage to provide the heating voltage for electron tubes. Continuously adjustable high-voltage source, which is safe to touch, with passive current limiting and digital voltage display.

- Current limiting to protect against contact with high voltages
- Safety transformer conforming to EN 61558-2-6
- Safe isolation between power supply and output circuits

High-voltage output: 0 – 5000 V DC, max. 2 mA, floating
 Heater voltage output: 6.3 V AC, max. 3 A, resistance to voltage up to 5 kV

1013412 648,90€



1008535

POWER SUPPLIES

➤ WITH 3B PRODUCTS - SAFE EXPERIMENTING IN CLASS!



Transformer 12 V, 25 VA (230 V, 50/60 Hz)

Simple transformer for student exercises. Short circuit proof, with connection leads and two cascadable 4 mm safety plugs.

- Safety extra-low voltage (SELV) and functional extra-low voltage (FELV)
- Safety transformer conforming to EN 61558-2-6
- Safe isolation between power supply and output circuits

Output: 12 V AC, max. 2 A
 Dimensions: 110x95x65 mm³ approx.
1000866 40,32€



Especially safe!

Transformer 12 V, 60 VA (230 V, 50/60 Hz)

Powerful electronic power supply, e.g. for the operation of lamps as used in optics. Short circuit proof, with connection leads and two cascadable 4 mm safety plugs.

- Safety extra-low voltage (SELV) and functional extra-low voltage (FELV)
- Safety transformer conforming to EN 61558-2-6
- Safe isolation between power supply and output circuits

Output: 12 V AC, max. 5 A
 Dimensions: 140x90x60 mm³ approx.
1020595 34,20€



Ideal for student experiments!



1003316

Transformer with Rectifier 3/ 6/ 9/ 12 V, 3 A (230 V, 50/60 Hz)

Extra low voltage power supply with overload protection contained in plastic housing. Output voltage switchable in four stages.

- Safety extra-low voltage (SELV) and functional extra-low voltage (FELV)
- Safety transformer conforming to EN 61558-2-6
- Safe isolation between power supply and output circuits

AC output: 3/ 6/ 9/ 12 V, max. 3 A
 DC output: 3/ 6/ 9/ 12 V, max. 3 A
 Dimensions: 210x170x90 mm³ approx.
 Weight: 2.6 kg approx.
1003316 117,90€

Especially safe!

Transformer with Rectifier 2/ 4/ 6/ 8/ 10/ 12/ 14 V, 5 A (230 V, 50/60 Hz)

Safety isolating transformer with safety cut out contained in metal housing. Output voltage switchable in 7 steps

- Safety extra-low voltage (SELV) and functional extra-low voltage (FELV)
- Safety transformer conforming to EN 61558-2-6
- Safe isolation between power supply and output circuits

AC output: 2/ 4/ 6/ 8/ 10/ 12/ 14 V, max. 5 A
 DC output: 2/ 4/ 6/ 8/ 10/ 12/ 14 V, max. 5 A
 Dimensions: 260x140x130 mm³ approx.
 Weight: 3.1 kg approx.

1003558 235,80€



1003558

DC Power Supply 1.5 – 15 V, 1.5 A (230 V, 50/60 Hz)

Handy DC power supply contained in a sturdy metal housing. The output voltage is continuously adjustable and is displayed via an analogue measuring instrument. The output is short circuit proof and floating.

- Safety extra-low voltage (SELV) and functional extra-low voltage (FELV)
- Safety transformer conforming to EN 61558-2-6
- Safe isolation between power supply and output circuits

DC output: 1.5 – 15 V, max. 1.5 A
 Residual ripple: 10 mV
 Dimensions: 100x80x150 mm³ approx.
 Weight: 2 kg approx.

1003560 128,70€

Ideal for student experiments!



1003560

More Power Supplies at 3bscientific.com!



1003562

AC/DC Power Supply 0 – 20 V, 5 A (230 V, 50/60 Hz)

Power supply with adjustable and stabilised DC voltage and analogue voltage and current display for DC voltage. The DC voltage component features an automatically alternating voltage and current control and is protected against continuous short circuits. The AC voltage can be selected in eight steps, the output is protected by an overcurrent circuit breaker.

- Safety extra-low voltage (SELV) and functional extra-low voltage (FELV)
- Safety transformer conforming to EN 61558-2-6
- Safe isolation between power supply and output circuits

DC output: 0 – 20 V, 0 – 5 A
 AC output: 2, 4, 6, 8, 10, 12, 15, 20 V, max. 5 A
 Residual ripple: <10 mV
 Dimensions: 235x175x245 mm³ approx.
 Weight: 8 kg approx.

1003562 430,20€

MULTIMETERS

Analogue Multimeters ESCOLA

Clear moving-coil instrument in shock-resistant plastic casing with two mirrored linear scales and clearly distinguishable measuring ranges. Includes battery test function and display of charge status as well as electronic calibration of zero point to the centre of the scale for all DC current and voltage ranges.

Use of a measurement amplifier ensures the measured values are linear even for AC voltages of up to 40 kHz. Only an inexpensive 1.5 V battery element is needed for operation. Nevertheless the meter will work for several years after any change of battery with normal usage, since the current discharge when in operation is no more than 2.5 mA maximum.

Scale length:	80 mm
Operating voltage:	1 – 3.5 V DC
Battery type:	Mignon, AA, R6
Accuracy:	Class 2 (DC), class 3 (AC)
Dimensions:	100x150x50 mm ³ approx.
Weight:	300 g approx.



1013526

No need to change fuses

Analogue Multimeter ESCOLA 30

Permanently short-circuit-proof student measuring instrument for measuring voltage and current in the safety extra-low voltage range. The electronic overload protection is achieved without the use of an equipment fuse, therefore obviating any need to change fuses or order spares. The protective system nevertheless operates without any auxiliary energy and is guaranteed even when the battery is flat or no battery is present.

Direct and alternating voltage:	0.3 – 30 V, 5 ranges each
Direct and alternating current:	1 – 3000 mA, 5 ranges each
Instrument category:	CAT I, 30 V

1013526 122,40€



THE IDEAL METER FOR STUDENT EXPERIMENTS

- + Unmistakeable measurement readings
- + Only an inexpensive 1.5 V battery element is needed for operation
- + Full functionality guaranteed even when the battery is no longer fully charged
- + Lithium batteries with higher open-circuit voltage can also be used
- + Battery protected by automatic cut-off after approximately 50 mins.
- + Distinct difference between 0 V display and the equipment being switched off due to inherently different position of needle

CAT III
600 V



1013527

Analogue Multimeter ESCOLA 100

Meter for classroom and practical experiments to measure voltage and current up to 600 V or 10 A respectively. Also features audible continuity testing. Includes a fuse to guarantee safety up to CAT III. The separate terminal sockets for current and voltage permit connection of the instrument that allows for current as well as voltage to be measured without having to reconnect the measuring leads. When switching from one measuring range to another, the circuit is never broken. All current measuring ranges are overload-proof for long-term current of up to 10 A.

Direct and alternating voltage:	0.1 – 600 V, 9 ranges each
Direct and alternating current:	0.1 mA – 3000 mA, 11 ranges each
Instrument category:	CAT III, 600 V

1013527 150,30€



Electrical safety of measuring instruments for current and voltage are assessed according to measurement categories stipulated in IEC 61010-1:

CAT I or unstipulated: Approved for measurements in circuits which are not directly connected to the low voltage mains grid (e.g. batteries).

CAT II: Approved for measurements in circuits which are directly connected, by a mains lead and plug for instance, to the low voltage mains grid (e.g. household or office appliance and lab equipment).

CAT III: Approved for measurements in circuits which are part of a building's wiring installation (e.g. stationary consumers, distribution terminals, appliances connected directly to the distribution box).

CAT IV: Approved for measurements in circuits which are directly connected to the source of the low voltage mains (e.g. electricity meters, main service feed, primary excess voltage protection). Note: the closer measurement is to be made to the low-voltage mains installation, the higher the measuring category needs to be.

> NEW!



1017895

Demo Multimeter

Electronic meter featuring a double scale for analogue measurement of current and voltage in demonstration experiments. It can handle measurements of current and voltage values and also allows the zero point to be set up in the centre of the scale for measurement of DC quantities. Switching between measuring ranges does not break any circuits connected to the equipment. This means it is possible to carry out measurements on voltage converters, for example, without causing induction surges. Resistance R, conductance G, impedance Z and admittance Y can easily be determined as quotients of current and voltage measurements thanks to the non-interrupting switch capability without the need to change the wiring.

This equipment is protected by fuses and authorised for making measurements in circuits directly connected to the low-voltage mains via plugs (CAT II), i.e. for measurements on house-hold appliances, for example. The current measuring ranges are resistant to long-term overloading up to 10 A. The meter is suitable for use as a free-standing instrument or for setting up in training panel frames.

Voltage ranges: 0.1 - 600 V AC/DC, 9 ranges
 Current ranges: 0.1 mA - 10 A AC/DC, 11 ranges
 Measuring category: CAT II: 600 V
 Dimensions: 259x297x125 mm³ approx.
 Weight: 1.7 kg approx.

1017895 536,40€



1002783

Digital Mini Multimeter

Very reasonably priced mini multimeter in pocket format for measuring voltage, DC current, resistance and temperature and also including diode and continuity tests. Overload protection for mA ranges, 10 amp range is unprotected. Includes measuring leads, type K thermocouple and battery.

DC voltage: 200 mV - 250 V, 5 ranges, $\pm 0.8\% \pm 2$ digits
 AC voltage: 200/ 250 V, 2 ranges, $\pm 1.2\% \pm 10$ digits
 DC current: 200 μ A - 10 A, 5 ranges, $\pm 1.0\% \pm 2$ digits
 Resistance: 200 Ω - 2000 k Ω , 5 ranges, $\pm 0.8\% \pm 2$ digits
 Temperature: 0 - 1000°C, $\pm 2.0\% \pm 3$ digits
 Safety classification: CAT II 250 V (IEC-1010-1)

1002783 14,40€

Digital Multimeter P1035

Compact 3½ digit multimeter for measuring voltage, current and resistance and also including diode and continuity tests. Complete with pouch, leads and battery.

DC voltage: 200 mV - 600 V, 5 ranges, $\pm 0.5\% \pm 2$ digits
 AC voltage: 200/ 600 V, 2 ranges, $\pm 1.2\% \pm 10$ digits
 DC current: 2000 μ A - 10 A, 4 ranges, $\pm 1\% \pm 2$ digits
 Resistance: 200 Ω - 2000 k Ω , 5 ranges, $\pm 0.8\% \pm 2$ digits
 Safety classification: CAT III 600 V (IEC-1010-1)

1002781 25,20€

Easy to replace fuses!



1002781



More Multimeters
at 3bscientific.com!

*10% discount for all products in this brochure by using action code, promotion closes 31/07/17

FUNCTION GENERATORS

+

ADVANTAGES

- + Easy and accurate adjustment
- + With built-in continuous sweep-mode
- + Ideal for recording resonance curves



1009957

Function Generator FG 100 (230 V, 50/60 Hz)

Function generator with power amplifier for use in versatile student and practical experiments covering simple harmonic oscillation, AC electricity and induction. Featuring illuminated, digital display for frequency, signal form, offset and other parameters. The output is short-circuit protected as well as being protected against induced voltages and spark discharges, e.g. for when experiment leads are unintentionally pulled out while coils are connected. In internal sweep mode, one trigger pulse is output per cycle and the voltage output is proportional to the frequency. With retractable feet. Includes power supply.

- Safety transformer conforming to EN 61558-2-6
- Safe isolation between power supply and output circuits

Frequency range: 0.001 Hz to 100 kHz
 Signal forms: Sine, square, triangular
 Offset: 0 to ± 5 V, adjustable in 0.1 V steps
 Output amplitude: 0 to 10 V, continuously adjustable
 Power output: 10 W, permanent
 Output current: 1 A, permanent, 2 A max.
 Sweep modes: External, continuous internal, individual internal
 Frequency range: 1 Hz to 100 kHz
 Time range: 0.04 s to 1000 s
 External sweep: Start via trigger pulse or application of 0 to 5 V control voltage
 Internal sweep: Start and stop via Start/Stop button
 One trigger output per cycle plus proportional voltage
 Power supply: Plug-in power supply, 12 V AC, 2 A
 Dimensions: 170x105x40 mm³ approx.
 Additional features: Fold-out feet

1009957 445,50€

+

ADVANTAGES

- + Ideal for introductory student experiments on AC
- + Simple generator for oscillations and waves



1017337

Function Generator SG10 (230 V, 50/60 Hz)

Sine-wave generator, which is particularly easy to use, featuring a power amplifier for use in student experiments. Includes 12 V AC plug-in power supply. One red and one green LED indicate the positive and negative half-waves of the output voltage. Their brightness corresponds to the configured amplitude. The way the output signal changes over time can be traced by means of an analogue voltmeter with zero-point in the centre or by means of an oscilloscope. The output is protected against short-circuits and against induced voltages as well as spark discharges.

- Safety transformer conforming to EN 61558-2-6
- Safe isolation between power supply and output circuits

Signal form: Sine-wave
 Frequency range: 0.01 – 10 Hz
 Output amplitude: 1 – 10 Vpp, continuously adjustable
 LED display: As of 2 V output voltage
 Output power: 1.5 W permanent
 Output current: 300 mA max.
 Distortion factor: <5%
 Connectors: 4-mm safety sockets
 Power supply: 12 V AC, 500 mA plug-in power supply

1017337 90,00€



Introductory experiment on AC

> NEW!



1020910

Digital Oscilloscope 2x30 MHz

Latest generation, dual-channel, digital storage oscilloscope with high-resolution colour display and large internal data memory.

- Mathematical functions including fast Fourier transforms (FFT)
- 20 Automatic measuring modes
- User-friendly operation featuring autoset and autoscale
- PASS/FAIL function implemented
- VGA output for connection to an external monitor
- LAN connection for remote connection via network
- USB connection for real-time data transmission or reading of internal memory

Includes two probes, two BNC cables, USB connecting cable, and software CD for Windows 2000/XP/VISTA/7/8/10.

1020910 413,25€

> NEW!



1020857

PC Oscilloscope, 2x25 MHz

Dual-channel, PC-based oscilloscope to be connected to a computer. Features powerful PC software for control and data requests. The highest possible degree of safety for users and the computer system thanks to galvanic isolation of the USB port.

- Mathematical functions including fast Fourier transforms (FFT)
- 20 Automatic measuring modes
- Data export for further processing (bin, txt, csv or xls)
- Image file for screenshots (png, bmp or gif)

Includes two probes (1:1, 10:1), a USB connecting cable, two BNC cables, operating instructions and software CD for Windows® XP/Vista/7/8

1020857 202,50€

	Digital Oscilloscope (1020910)	PC Oscilloscope (1020857)
Channels:		Two
Band width:	30 MHz	25 MHz
Sample rate:	250 MSa/s	200 MSa/s
Operating modes:		CH1, CH2, XY
Input coupling:		DC, AC, GND
Input impedance:		1 MΩ ±2% 10 pF ± 5 pF
Input voltage:		0 – 400 V DC or ACpp
Deflection coefficient:	2 mV/div. – 10 V/div.	2 mV/div. – 50 V/div.
Time-base coefficient:		5 ns/div. – 100 s/div.
Trigger type:	Edge, Video, Pulse, Slope	Alternate, Edge, Video, Pulse, Slope
Trigger modes:		Auto, Normal, Single
Trigger detection:		Sample, Peak Detect, Average
Memory size:	10000 measurements	5000 measurements
Interface:	USB 2.0, VGA, LAN	USB 2.0
Power supply:	100 – 240 V, 50/60 Hz	via two USB ports
Dimensions:	355x178x118 mm ³ approx.	170x120x18 mm ³ approx.
Weight:	1,6 kg approx.	260 g approx.

DIGITAL THERMOMETERS

> NEW!



1020908

Thermal Imaging Camera

Modern infra-red thermal imaging camera for producing images of infra-red radiation from an object based on detected infra-red radiation in relation to the ambient temperature.

- User-friendly graphic menu operation
- Photography using built-in digital camera
- Up to 25000 photos can be saved on Micro SD card
- Recordings featuring time and date documentation
- Images with emission factor and measurements
- Five colour palettes for thermal imaging
- Five levels of photograph and thermal imaging superimposition
- Cross-hairs, plus cold-spot and hot-spot display
- Minimum and maximum value display
- Automatic shut-off

Includes case, batteries, Micro SD card and instruction manual.

Temperature range:	-20°C ... 300°C
Resolution:	0.1°
Sensitivity:	0.3°C
Precision:	±2% or 2°C
Display:	60 mm (2.4") LCD-TFT
Thermal image resolution:	60 x 60 pixels
Field of vision:	20° x 20°
Emission factor:	Adjustable from 0.1 – 1.0
Wavelength:	8 – 14 μm
Image frequency:	6 Hz
Focus range:	50 cm (fixed)
Memory:	Micro SD card
Voltage supply:	4 x 1.5 V AA batteries
Display:	Multi-line, multi-function display
Dimensions:	212 x 95 x 62 mm ³ approx.
Weight:	320 g approx.

1020908 474,05€

> NEW!



1002791



1020909



ADVANTAGES

- + Fast, easy and accurate measurements
- + Automatic selection of measurement range
- + Practical single-handed operation
- + Modern, handy design
- + Large-scale 3½-digit LCD display

Infra-red Thermometers

Surface thermometer for contactless temperature measurement from a safe distance, e.g. in inaccessible places, hot or moving objects. With laser diode for laser sighting, illuminated LCD display, range overflow display, measured value storage function, selection between Celsius and Fahrenheit, automatic switch off.

The infra-red thermometer 380°C D permits rapid measurement of temperature differential with the LED display (red, green or blue). Including case, battery and instruction manual.

	1002791	1020909
Designation	Infra-red Thermometer, 800°C	Infra-red Thermometer, 380°C D
Measuring range	-50° C – +800° C	-50° C – +380° C
Accuracy	±1 % of measured value	±2 % of measured value
Response time	150 ms	< 1 s
Optical resolution	20:1	10:1
Max. temperature display	yes	---
Voltage supply	9 V battery	9 V battery
Dimensions	146x43x104 mm ³ approx.	200x124x50 mm ³ approx.
Weight	170 g approx.	220 g approx.
Price	115,20€	41,33€

Prices are valid until 31/12/17, VAT and postage costs not included.

Digital Thermometer Type K/IR

Digital two channel thermometer with two K-type inputs and additional external infra-red sensor. Can also be used for measurements at low temperatures. With automatic shut off, maximum value storage and data hold function. Includes case, 2 K-type thermocouple sensors, infra-red temperature sensor, 9 V battery and instruction manual.

Measurement inputs: 2x K-type, external IR input
 Measuring functions: T1, T2, T3, T1-T2, T1-T3, T2-T3
 Measuring range: -200 – 1372°C (type K), -30 – 550°C (IR)
 Measurement error: ±0.5% + 2°C (type K), ±2.5% + 2°C (IR)
 Resolution: 0.1°C
 Unit of measurement: °C or K
 Emission factor: 0.95 fixed
 Digital display: 3¾ digit LCD
 Background lighting: blue
 Voltage supply: 9 V battery
 Dimensions: 75x200x50 mm³ approx.
 Mass: 280 g approx.

1002799 151,20€



1002799

Infrared Temperature and Humidity Gauge

Digital measuring device for contact-free temperature measurement from large distances, e.g. of hot or moving objects or inaccessible points of measurement, and for simultaneous humidity display. With laser diode as detection aid, integrated in the measuring probe, illuminated LCD display, max and data-hold function, switchable between °C and °F, automatic switch-off. Includes pouch and battery.

Measuring range, temperature: -50° C to +500° C
 Divisions: 0.1° C
 Accuracy: ± 2% of measured value ± 2° C
 Measuring range, humidity: 5% to 95%
 Divisions: 0.1%
 Accuracy: ± 3.5%
 LCD dual-function display: 3 ½-digit, 21 mm with backlighting
 Voltage supply: 9 V battery
 Dimensions: 90x170x45 mm³ approx.
 Mass: 360 g approx.

1002795 82,80€



1002795

Digital Thermometers

Versatile digital thermometers for type-K temperature sensors with single or dual input (1002794) for measuring instantaneous or differential temperature (T1 – T2 1002794). With storage of maxima and Data-Hold function. Includes type-K temperature sensor (1002794 2x), battery, holster and carrying bag.



	1002793	1002794
Designation	Digital Thermometer, 1 Channel	Digital Thermometer, 2 Channels
Measuring range	-50° C – +1300° C 223 K – 2000 K	-50°C – +1300°C
Division	0.1° C, 1 K	0.1°C/F
Accuracy	±0.5% +1° C	±0.5% +1°
Display	3½ digit illuminated LCD	3½ digit illuminated LCD
Digit size	21 mm	21 mm
Voltage supply	9 V battery	9 V battery
Dimensions	90x170x45 mm ³ approx.	90x170x45 mm ³ approx.
Weight	350 g approx.	350 g approx.
Price	52,20€	69,30€

SOUND LEVEL METERS



1002778



1002780

Sound Level Meter P5055

Digital measuring instrument for universal application in detecting noise levels from a variety of sound sources over a broad range. Features built-in calibration signal plus maximum value and value hold functions. Slow mode for average noise level and fast mode for recording brief sound sequences and determining maximum noise level. A-weighting of frequency (based on human hearing) for open-air measurements and also C-weighting, e.g. for measurements of engine noise. Robust plastic casing, analogue output for external measuring instruments, threaded hole for mounting on a stand. Foam-filled carry case.

Measurement range: 35 – 130 dB
 Resolution: 0.1 dB
 Accuracy: ±3.5 dB at 94 dB (1kHz)
 Display: 3½-digit LCD, 17 mm
 Microphone: Electret capacitor microphone
 Voltage supply: 9 V block-type battery
 Dimensions: 251x64x40 mm³ approx.
 Weight: 250 g approx.

1002778 84,60€



1012741

Noise Level Meter P8005

Digital noise meter with background noise suppression for all types of measurements of ambient noise, e.g. for measuring noise levels in schools, offices, factories, traffic and homes or for noise projects. Includes data logger and USB port for long-term measurements. Choice of manual and automatic operating modes. Capability for min. and max. measurements.

Includes case, USB cable, Windows software, stand, 9 V mains adaptor, 9 V battery and instruction manual.

Frequency range: 31.5 Hz – 8 kHz
 Dynamic range: 50 dB
 Level ranges: 30 – 80 dB (low)
 50 – 100 dB (medium)
 80 – 130 dB (high)
 30 – 130 dB (automatic)

Precision: ±1.4 dB
 Digital display: 4 digit LCD, 20 mm
 Multi-functions display: Digital display of measurement, measuring time, bar graphs plus overs and unders

Applicable standards: IEC-61672-1 type 2, ANSI S1.4 type 2
 Response times: 125 ms (fast), 1s (slow)
 Microphone: ½-inch, with electret capacitor
 Analogue output: AC/DC
 Voltage supply: 9 V battery or 9 V mains adaptor
 Dimensions: 90x280x50 mm³ approx.
 Weight: 350 g approx.

1002780 208,80€



NOISE ...

- ... damages hearing
- ... makes it more difficult to hear genuinely important signals
- ... impedes both physical and mental work
- ... disturbs and adversely affects well being
- ... disturbs relaxation and sleep
- ... can cause chronic stress, physical ailments and illness

Noise Level Indicator SPL

Handy and easy-to-use noise level meter with digital display in decibels (dB) and an arbitrarily adjustable trigger threshold for use as a traffic-light style noise indicator with a happy green face and a sad red face. Can be mounted on a wall or set up on a table top. Its well-conceived compact design makes it easy to transport. Automatically switches to electricity-saving stand-by mode when noise is low for a long period. The brightness of the display can also be adjusted. Includes a stand base, USB/miniUSB cable and USB power supply.

Display: 100 mm diam, with LED
 Measuring range: 40 dB to 130 dB
 Resolution: 1dB
 Thresholds for colour display: Adjustable to any level in steps of 1 dB
 Voltage supply: 5 V DC via miniUSB socket
 Power consumption: 150 mA (normal operating mode)
 <1 mA (stand-by)

USB power supply: 100 – 240 V, 50/60 Hz
 Dimensions: 130x145x12 mm³ approx.
 Weight: 400 g approx.

1012741 302,40€

pH Meter (2 in 1)

Digital pH meter for the simultaneous measurement of the pH value of aqueous fluids and their temperatures. The measurement of the pH value is performed by determining the electrical potential difference between acidic, neutral and base fluids. Robust, water-proof housing, large LCD display with permanent background illumination and simple to operate. Including calibrating solution, batteries and instruction manual.

pH range:	0 – 14 pH
Resolution:	0.01 pH
Accuracy:	± 0.05 pH
Temperature compensation:	0 – 50°C
Temperature measurement:	0 – 50°C
Resolution:	0.1°C
Accuracy:	± 1°C
Voltage supply:	4 x 1.5 V (AG-13) batteries
Protection class:	IP 65 water-proof
Display:	3½-digit LCD display, 11 mm, max. 1999
Dimensions:	190x35x35 mm ³ approx.
Weight:	100 g approx.

1020914 55,10€



NEW!

1020914

1020915

pH Meter

Digital pH measuring instrument for the measurement of the pH value of aqueous liquids using the immersion probe tips to determine the electrical potential difference between acidic, neutral and basic liquids. The device has a robust housing with compact dimensions and is easy to operate. It is equipped with a large LCD display with continuous background illumination, 2 adjustment potentiometers for calibrating to pH = 4 or to pH = 7 using the matching screwdriver. Including calibration solution, screwdriver, battery and instruction manual.

pH range:	0 – 14 pH
Resolution:	0.01 pH
Accuracy:	± 0.05 pH
Temperature compensation:	0 – 50°C
Voltage supply:	9 V battery (NEDA 1604)
Display:	3½-digit LCD display, 18 mm, max. 1999
Dimensions:	150x70x25 mm ³ approx.
Weight:	230 g approx.

1020915 116,85€

NEW!



1020907

Laser Range Finder

Professional laser range finding instrument with multi-lined LCD display and background illumination especially designed for distance measurements of extremely high precision and for locations difficult to access. Speed buttons for direct and indirect measurement (according to Pythagoras), area and volume calculation, addition and subtraction operations. With internal memory for 99 recorded measurement values, retractable 90° bracket for precise targeting of the measurement point, spirit level and tripod socket. Including case, batteries and instruction manual.

Measurement range:	0.05 – 60 m
Measurement units:	m (metre), in (inch), ft (feet)
Accuracy:	± 2 mm
Internal memory:	99 values
Laser:	620 nm – 680 nm, <1 mW, class: 2
Voltage supply:	2 x 1.5 V AAA batteries
Display:	Multi-lined multifunction display
Dimensions:	118x54x28 mm ³ approx.
Weight:	135 g approx.

1020907 76,95€

Digital Luxmeter

Reasonably priced, easy to use pocket lux-meter for testing and measurement of light conditions. C.I.E. standard spectrum. Including light sensor, pouch and battery.

Measuring ranges:	200 – 50000 lux, 4 ranges, ±5%
Voltage supply:	12 V battery (A23)
Dimensions:	65x115x25 mm ³ approx.
Weight:	160 g approx.

1002779 35,10€



1002779

GEIGER COUNTERS



1002722



1012894

Geiger Counter

Versatile, easy to use and compact precision instrument for measuring α -, β - and γ -radiation. With filter selection switch for filtering out types of radiation, large display and integrated USB interface. Including USB cable, Windows software, and operating instructions.

The following functions and operating modes are available for measurement:

- Standard mode for displaying the current radiation level. Also equipped with variable acoustic and optical warning threshold signal and display of average radiation from previous day.
- Pulse counting either permanent or with variable gate time. Additional optional acoustic count indication.
- Count rate measurement.
- Integrated display of date and time.
- The number of pulses registered is stored in the internal memory. This facilitates recording e.g. of weekly values for up to 10 years.
- Computer docking station. The software enables the measured data to be evaluated and processed on an MS-Windows PC.

Radiation types:	α from 4 MeV, β from 0.2 MeV, γ from 0.02 MeV
Measured variables:	equivalent dose in Sv/h, mSv/h, μ Sv/h pulses/s, pulses/variable time interval
Display:	LCD, 4 digit, numerical with display of measured variable, quasi analogue bar chart, operating mode indicators
Radiation detector:	End window Geiger-Müller counter tube, stainless steel housing with neon-halogen filling
Measuring length:	38.1 mm
Measuring diameter:	9.1 mm
Mica window:	1.5 – 2 mg/cm ²
Gamma sensitivity:	114 pulses/min for ⁶⁰ Co radiation = 1 μ Sv/h in background radiation energy band
Background rate:	10 pulses per minute approx.
Internal memory:	2 kilobytes
Battery life:	3 years approx.
	1002722 344,70€

Dosimeter Radex RD 1706

Used for determining dose rates in μ Sv/h for β -, γ - and X-rays, this radiation meter can be operated by non-professionals while nonetheless offering the features of a professional dosimeter. Including two built-in Geiger-Müller counter tubes and a large, illuminated LCD display. The device measures the activity of β - and γ -particles and uses the results to calculate the dose rate. Detection of each particle is indicated by an audio signal to facilitate searching for radioactive sources. The difference between the mean dose rate and background radiation level, as well as the background radiation level itself are displayed in the „background“ mode. Measured values remain saved after the device has been turned off.

Counters:	Two GM counter tubes SBM20-1
Measurement variable:	Ambient equivalent dose rate H*(10)
Measuring range:	0.05 ... 999.0 μ Sv/h
Alarm threshold:	Adjustable from 0.10 to 99.0 μ Sv/h
Alarm:	Audio or vibration signal
Measurement and calculation times:	26 s, 1 s (at H*(10) > 3.5 μ Sv/h)
Value display duration:	Continuous
Energy detection range	
X-radiation and γ -radiation:	0.03 to 3.0 MeV
β -radiation:	0.25 to 3.5 MeV
Batteries:	1.5 V, AAA (1 x or 2 x)
Operating time:	500 h
	1012894 250,20€

Geiger-Müller Counter Tube

Self-quenching halogen pulse ionisation chamber for detecting alpha, beta, gamma and x-ray radiation. In metal housing with mica window, removable mounting clamp with shaft. Long plateau length.

Filling:	Neon/argon mixture, halogen as quenching agent
Cathode dimensions:	approx. 39x14 mm ²
Window:	mica, 9 mm dia.
Mass per unit area:	1.5 – 2.0 mg/cm ²
Plateau length:	400 V – 600 V
Operating voltage:	400 – 600 V (recommended: 500 V)
Relative plateau slope:	0.04 %/V
Dead time:	90 μ s
Limiting resistor:	10 M Ω , integrated in holder
Shaft:	100 mm x 10 mm dia. approx.
	1001035 213,30€

Additionally required:
1002746 HF Patch Cord, 1 m 15,84€
1001033 Digital Counter (230 V, 50/60 Hz) 347,40€




1001035

1001033



> NEW!



Electrometer (230 V, 50/60 Hz)

Impedance converter with high-resistance input for measuring extremely small charges and currents. The input signal is converted into a proportional voltage, which can then be measured with an external voltmeter. During the measurement the potentials of the electrometer and the experimenter must be matched by using a metal rod connected to earth. Includes a 12 V AC plug-in power supply.

Electrometer gain factor:	1.00
Input resistance:	$>10^{12} \Omega$
Output resistance:	$<1 \text{ k}\Omega$
Input current:	$<10 \text{ pA}$
Input capacitance:	$<50 \text{ pF}$
Max. output voltage:	$\pm 10 \text{ V}$
Resistance to excess voltage:	1 kV (from low-resistance sources) 10 kV (from high-resistance sources)
Supply voltage:	12 V AC
Dimensions:	110x170x30 mm ³ approx.
Weight:	1 kg approx.

1001025 157,50€

Additionally recommended:

- 1006813 Electrometer Accessories**
- 1013526 Analogue Multimeter Escola 30**
- 1008535 DC Power Supply 450 V (230 V, 50/60 Hz)**

Electrometer Accessories

Set of accessories for carrying out basic experiments on electrostatics, electricity and the photoelectric effect in combination with an electrometer (1001025) and 450 V DC power supply (1008535).

Contents:

- 1 Faraday cup
- 1 Pair of friction rods
- 1 Metal rod with 4-mm drilled hole
- 1 Safety adaptor socket
- 1 Plug-in capacitor 1 nF
- 1 Plug-in capacitor 10 nF
- 1 Plug-in resistor 100 M Ω
- 1 Plug-in resistor 1 G Ω
- 1 Plug-in resistor 10 G Ω
- 1 Zinc electrode
- 1 Grid electrode

1006813 171,90€



Measurement Amplifier U (230 V, 50/60 Hz)

Measurement amplifier U amplifies low-amplitude measurement signals from low-resistance signal sources for measurement with any chosen voltmeter or oscilloscope. By using an external shunt resistor it is also possible to measure small currents. Offset voltages can be compensated using coarse and fine offset adjustment knobs. Amplification (gain) can be selected in ranges from 0 to 5 powers of ten. High-frequency noise or other interference signals are filtered out by means of a low-pass filter with step-wise selectable time constants between 0 and 3 seconds. The output voltage has the same sign as the input voltage.

Input resistance:	10 k Ω
Output resistance:	300 Ω
Offset voltage drift:	$< 2 \mu\text{V/K}$ (after 15 mins. operation approx.)
Gain factors:	10 ⁰ ; 10 ¹ ; 10 ² ; 10 ³ ; 10 ⁴ ; 10 ⁵
Tolerance for gain factors:	$< 2.5 \%$
Input voltage:	max. $\pm 12 \text{ V}$ (overload protected for brief transients up to 100 V)
Output voltage:	0 ... $\pm 12 \text{ V}$ (short-circuit protected)
Power supply (via plug-in supply provided):	12 V AC
Ambient temperature:	5°C ... 23°C ... 40°C
Storage temperature:	-20 ... 70°C
Relative humidity:	$< 85\%$ no condensation
Operational alignment:	Horizontal
Contamination level:	2
Protection class:	IP20
Dimensions:	170x105x50 mm ³ approx.
Weight:	335 g approx.

1020742 315,90€



3B Scientific

A worldwide group of companies

UK 3B Scientific Limited

Suite 1, Formal House, Oldmixon Crescent


Weston-super-Mare • BS24 9AY

Phone: 01934 425333 • Fax: 01934 425334

3bscientific.com • uk3bs@3bscientific.com



Our physics experiments can be found online at 3bscientific.com or ask for our catalogues.



UE8020250: Investigation of an island grid or microgrid used to generate and store electrical energy

*10% discount for all products in this brochure by using action code, promotion closes 31/07/17