

Build your own Labor and Delivery Suite.



Operative Experience

RealMom	
Vaginal Delivery 1021404	<ul style="list-style-type: none"> + Full body female simulator + Full-term baby + Instructor Tablet + Patient Display Monitor
C-Celia	
Fetal Extraction 1021386	<ul style="list-style-type: none"> + Open abdomen simulator and uterus + 1 large baby and 1 set of twin babies + Supports thousands of deliveries
C-Section Skills 1021392	<ul style="list-style-type: none"> + Abdominal and uterine incision repair simulator + Pregnant abdomen and uterus + Flaccid abdomen and empty uterus + 50 abdominal and 50 uterine inserts with amniotic sacs
C-Section Delivery 1021391	<ul style="list-style-type: none"> + Pregnant abdomen simulator with pressurized bleeding & lubrication system + 5 loadable uteri + 13 abdominal inserts + 1 large baby and 1 set of twin babies
Postpartum Hemorrhage 1021387	<ul style="list-style-type: none"> + Open abdomen simulator with pressurized bleeding & lubrication + 2 multi-use uteri
Emergency Hysterectomy 1021389	<ul style="list-style-type: none"> + Open abdomen simulator with pressurized bleeding & lubrication + Uterus, bladder, ureters and detailed representations of uterine arteries

Available Replacement Parts


- Uterine incision inserts for C-Section trainer
1021393
- Abdominal inserts and repair inserts
(each enables two incisions and two repairs)
1021394
- Emergency hysterectomy replaceable uteri
1021390
- Replaceable abdominal insert and pregnant uterus
1021396

LABOR AND DELIVERY SUITE



The World's Only Complete Labor and Delivery Suite

Introducing the first-ever, integrated solution for labor and delivery training, featuring a **complete suite of fully operable, hands-in-the-body simulators** for delivery, C-Section and emergency obstetric procedures.



3B Scientific GmbH • Rudorffweg 8 • 21031 Hamburg • Germany • 3bscientific.com
Phone: +49 40 73966-0 • Fax: +49 40 73966-100 • E-mail: info@3bscientific.com

3bscientific.com

The Labor and Delivery Suite from Operative Experience.

Finally, a complete and integrated training solution for the entire labor and delivery continuum of care.

Say goodbye to the old, unreliable plastic manikin technology and discover what real, hands-in-the-body simulation looks like. OEI simulators are fully operable and feature incredibly accurate and tactile anatomy, realistic blood and fluids, and are uniquely designed to provide learners with a high-adrenaline, real-world operative experience.

And, with OEI, there's no complicated programming or setups and no simulator down time – because our products simply work.

RealMom™



Vaginal Delivery

The world's most natural and realistic vaginal delivery simulator.

- + Natural delivery with active, soft tissue birth canal
- + Realistic, soft tissue and accurate internal anatomy
- + Full-body female simulator and full-term baby with unsurpassed realism and lifelike appearance
- + Provides a full-range of delivery procedures including forceps and vacuum assist
- + Postpartum hemorrhage control
- + Instructor control of dilation and delivery progression
- + Physiological presentation and control of maternal and fetal heart rate, fetal heart tones and contractions, pulses, blood pressure and SPO₂

C-Celia™



C-Section

A complete C-Section delivery training system available as a full suite or in custom configurations designed to meet your learning needs.

- + Fetal Extraction Simulator – open, post C-Section abdomen enabling difficult delivery maneuvers including forceps and suction assist
- + C-Section Skills Trainer – abdominal and incision repair simulator with pregnant abdomen and uterus for incision practice and tissue repair
- + Emergency C-Section Delivery Simulator – pregnant abdomen simulator enabling full transverse or vertical emergency C-Section through abdominal wall and uterus with baby, placenta, and pressurized lubrication system with simulated amniotic fluid



Postpartum Hemorrhage

Open abdomen simulator for practicing hemorrhage control techniques.

- + Postpartum hemorrhage control trainer with integrated pressurized bleeding and lubrication for repair of complex lacerations and uterine artery injuries
- + Open abdomen with empty, durable, multi-use uterus enables scores of procedures before uterus placement
- + Supports hemorrhage control techniques including Bakri Balloon, B-Lynch suture, Hayman suture and O'Leary Stitch



Emergency Hysterectomy

Training for postpartum and caesarean hysterectomy procedures including hemorrhage control and fully-operable removal of the uterus.

- + Open abdomen with bleeding and integrated lubrication includes durable, flaccid, abdominal wall insert, attachable, partially contracted uterus
- + Model contains a uterus, bladder, ureters and detailed representations of uterine arteries