

1020592 | P100

3bscientific.com

Trainer for wound care and bandaging techniques

TRAINER FOR WOUND CARE AND BANDAGING TECHNIQUES



Introduction


You have purchased a wound care and bandaging technique trainer with which you can demonstrate professional wound care, and on which students can practice. The trainer has 13 wound sites that can be discussed and cared for. Standards (e.g. hygiene) as part of professional wound management in nurse and medicine can be practiced and improved.

Delivery contents

Torso with 13 wound sites
Carry bag

Functions

The trainer for wound care and bandaging techniques is for the following surgical procedures and physical injuries:



The wounds and drain systems shown can be assessed professionally and cared for. Wound care can be carried out using the usual means and therefore provides training situations that are as close to reality as possible. Given the flexible structure of the material used, students can practice putting on and changing bandages and dressings in realistic conditions.

Cleaning and care

The surfaces of the simulator can be cleaned with a damp cloth. Any residue from the plaster bandages can be removed with normal plaster remover. Please ensure that the trainer is completely dry before storing.

Technical data

Dimensions:
Simulator (H x W x D): 23 x 46 x 77 cm
Weight: approx. 4.3 kg
Operating conditions
Operating temperature: -10°C to +40°C
Storage temperature: -20°C to +40°C



3B Scientific

A worldwide group of companies

3B Scientific GmbH • Customer Service Medical
Rudorffweg 8 • 21031 Hamburg • Germany
Phone: + 49 (0)40-73966-0 • Fax: + 49 (0)40-73966-100
3bscientific.com • info@3bscientific.com

© Copyright 2016 for instruction manual and design of
product: 3B Scientific GmbH, Germany