

STORABLE MAT

low height, high performance



specifications

DYNAMIC LOAD	550 lbs.
WEIGHT	40-45 lbs
HEIGHT RANGE	16" - 25"
WIDTH	33", 35", 40"
LENGTH	78"
UPHOLSTERY	TerraTouch™ , Tufflex™
TOP	Rectangular
PADDING	Semi-Firm , Firm
<i>Some configurations may not be available</i>	
Bold - Standard	<i>Italic - Additional Charge</i>

The only portable mat table on the market with a 550 lb weight capacity, the OAKWORKS® portable mat is as flexible as it is stable. Perfect for space constrained facilities where mat work is required, the portable mat table provides up to a 40" by 78" working surface while easily stores out of the way in a super-compact suitcase sized space.

All Portable Tables :U.S. Patent #5,009,170
 CS System™ (Complementary Suspension) : U.S. Patent #6,192,809
 UniLock™: U.S. Patent # 7,472,656



**BARIATRIC RATED
 550 LB DYNAMIC
 LOAD CAPACITY**

**ULTRA STIFF LEG BRACING
 HOLDS UP UNDER HEAVY
 LATERAL FORCE**

Details Make the Difference



1.

2.

3.

1. INTEGRA-HINGE™

Ingenious patented design reinforces the table to support heavy loading which no piano hinge can match.

2. STABILIZER FOOT

For maximum stability each 40" Storable Mat comes standard with stabilizing feet, which easily fold out of the way for storage.

3. QUICK PORTABILITY

With its space saving design, the Storable Mat folds up easily for storage and is ready to follow wherever you go.



Favorite Manufacturer Support - 2009
 Favorite Treatment Table Manufacturer - '07, '08, '09



800.558.8850

www.massagetables.com
 MMCSL0002 rev.42010

UPHOLSTERY

TERRATOUGH™

One more investment in wellbeing.

For years, vinyl, also known as polyvinyl chloride, has been the standard for massage tables. OAKWORKS® is proud to bring you our TerraTouch™, a revolutionary, earth-friendly upholstery fabric, made from environmentally friendly, PVC-free premium polyurethane. We designed it with the best features for bodywork.

- a soft “hand” for ultra comfort
- a great memory to retain its shape
- excellent abrasion and tear resistance
- high oil and cold endurance
- easy to clean



TUFFLEX™

This durable fabric is our abrasion-resistant upholstery, ideal for Athletic Trainers and Sports Therapists.



Favorite Manufacturer Support - 2009
Favorite Treatment Table Manufacturer - '07, '08, '09



PADDING

AERO • CEL™ PADDING

New technology in padding options now gives you 4” of more resilient support with a softness that conforms to the body’s curves. Offering the most relaxing and nurturing environment possible, it’s an investment that will return long years of “ahhhhs”.

Weighs 4-5 pounds more than semi-firm.



PLUSH PADDING

Our traditional 3” multi-layer, super supportive multi-density system for optimal comfort and durability. Its perfect balance of softness and support make this padding ideal for bringing relaxation to clients everywhere.

Weighs 2-3 pounds more than semi-firm.



SEMI-FIRM PADDING

2 1/4” multi-layer, multi-density system for firmer support and comfort with easier portability. Designed for the best balance of lighter weight and comfortable support. We use semi-firm on our portable tables to determine their published weights.



FIRM PADDING

1 1/2” multi-layer, multi-density system provides firm support. Designed for Physical Therapy and Athletic Training needs. Weighs 1-2 pounds more than semi-firm.



800.558.8850
www.massagetables.com